Datasheet revision 1.1 www.chipquik.com

Flux Remover (Non-Flammable) in 30ml (1.0oz) Squeeze Bottle w/tip

Product Highlights

A Flux Remover that is Non-Flammable
Can be used with Leaded and Lead-Free applications
RoHS 3 and REACH compliant

CHIPQUIK CONTRM-1.0 Sa Immover (Non-Farmis Control 1.0 Sa Immove

Specifications

Packaging: 30ml Squeeze Bottle

Shelf Life: Refrigerated >24 months, Unrefrigerated >24 months

Storage and Handling

Store refrigerated or at room temperature 3-25°C (37-77°F). The flux remover can be allowed to freeze. Freezing will not degrade this product. Allow 4 hours for flux remover to reach an operating temperature of 20-25°C (68-77°F) before use.

Transportation

This product has no shipping restrictions. Shipping below 0°C (32°F) or above 25°C (77°F) for normal transit times by ground or air will not impact this product's stated shelf life.

Conforms to the following Industry Standards:

RoHS 3 Directive (EU) 2015/863:

Yes

Flux Remover Selection Guide

Description	Packaging	Orderable Part Number
Degreaser / Contact Cleaner (Non-Flammable)	Pen, 10ml	CQNRM
	Squeeze Bottle, 0.5oz	CQNRM-0.5
	Squeeze Bottle, 1.0oz	CQNRM-1.0
Flux Remover (Non-Flammable)	Pen, 10ml	CQNTRM
	Squeeze Bottle, 0.5oz	CQNTRM-0.5
	Squeeze Bottle, 1.0oz	CQNTRM-1.0
Flux Remover (Fast Evaporating, Non-Flammable)	Pen, 10ml	CQNIRM
	Squeeze Bottle, 0.5oz	CQNIRM-0.5
	Squeeze Bottle, 1.0oz	CQNIRM-1.0
Flux Remover (Foam-Free, Flammable)	Pen, 10ml	CQIRM
	Squeeze Bottle, 0.5oz	CQIRM-0.5
	Squeeze Bottle, 1.0oz	CQIRM-1.0
Flux Remover (Foaming, Flammable)	Pen, 10ml	CQ6RM
	Squeeze Bottle, 0.5oz	CQ6RM-0.5
	Squeeze Bottle, 1.0oz	CQ6RM-1.0